

TBC1D4 antibody

Product Information

Catalog No.:	FNab09889
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


May act as a GTPase-activating protein for RAB2A, RAB8A, RAB10 and RAB14. Isoform 2 promotes insulin-induced glucose transporter SLC2A4/GLUT4 translocation at the plasma membrane, thus increasing glucose uptake.

Immunogen information

Immunogen:	AS160
Synonyms:	TBC1 domain family member 4 Akt substrate of 160 kDa (AS160) TBC1D4 AS160 KIAA0603
Observed MW:	147 kDa
Uniprot ID :	O60343

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	



RAW 264.7 were subjected to SDS PAGE followed by western blot with FNab09889(AS160 antibody) at dilution of 1:1000