
Goat anti- human IgA

Product Information

Catalog No.:	FNab09895
Size:	1mg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Goat
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	human IgA
Synonyms:	IgA1, IGHA1
Observed MW:	55-60 kDa

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:1000-1:10000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



human saliva tissue were subjected to SDS PAGE followed by western blot with FNab09895(Human IgA Antibody) at dilution of 1:5000