

## p24 antibody

### Product Information

Catalog No.:	FNab09919
Size:	100µg
Form:	liquid
Purification:	protein A+G purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	5C6
IsoType:	IgG1
Storage:	0.01MPBS pH 7.4, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	Recombinant Capsid protein p24
Synonyms:	Gag polyprotein Pr55Gag Matrix protein p17 (MA) Capsid protein p24 (CA) Spacer peptide 1 (SP1) Alternative names: p2 Nucleocapsid protein p7 (NC) Spacer peptide 2 (SP2) gag
Observed MW:	None
Uniprot ID :	P04591

### Application

Reactivity:	HIV-1
Tested Application:	ELISA