

CTLA4 antibody

Product Information

Catalog No.:	FNab09937
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	cytotoxic T-lymphocyte protein 4
Synonyms:	Cytotoxic T-lymphocyte protein 4 Cytotoxic T-lymphocyte-associated antigen 4 (CTLA-4) CTLA4 CD152
Observed MW:	37 kDa
Uniprot ID :	P16410

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:20 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human tonsil tissue slide using FNab09937 (CTLA4 antibody) at dilution of 1:100



mouse lung were subjected to SDS PAGE followed by western blot with FNab09937 (CTLA4 antibody) at dilution of 1:1000