

## AKT1S1 antibody

### Product Information

Catalog No.:	FNab09942
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	Proline-rich AKT1 substrate 1
Synonyms:	Proline-rich AKT1 substrate 1 40 kDa proline-rich AKT substrate AKT1S1 PRAS40
Observed MW:	50 kDa
Uniprot ID :	Q96B36

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:20 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab09942( AKT1S1 Antibody) at dilution of 1:100



mouse heart were subjected to SDS PAGE followed by western blot with FNab09942(AKT1S1 Antibody) at dilution of 1:1000