

RP2 antibody

Product Information

Catalog No.:	FNab09990
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Acts as a GTPase-activating protein(GAP) involved in trafficking between the Golgi and the ciliary membrane. Involved in localization of proteins, such as NPHP3, to the cilium membrane by inducing hydrolysis of GTP ARL3, leading to the release of UNC119(or UNC119B). Acts as a GTPase-activating protein(GAP) for tubulin in concert with tubulin-specific chaperone C, but does not enhance tubulin heterodimerization. Acts as guanine nucleotide dissociation inhibitor towards ADP-ribosylation factor-like proteins.

Immunogen information


Immunogen:	retinitis pigmentosa 2
Synonyms:	Protein XRP2 RP2
Observed MW:	40 kDa
Uniprot ID :	O75695

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human intestinal cancer using FNab09990(XRP2 antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab09990(XRP2 antibody) at dilution of 1:1000