

RHD antibody

Product Information

Catalog No.:	FNab10000
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	Blood group Rh(D) polypeptide
Synonyms:	Blood group Rh(D) polypeptide RHXIII Rh polypeptide 2 (RhPII) Rhesus D antigen RHD
Observed MW:	34 kDa/41 kDa/43 kDa/45 kDa/53 kDa
Uniprot ID :	Q02161

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human endometrium tissue using FNab10000(RHD Antibody) at dilution of 1:200



Rat heart were subjected to SDS PAGE followed by western blot with FNab10000(RHD antibody) at dilution of 1:1000