

EREG antibody

Product Information

Catalog No.:	FNab10012
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	Proepiregulin
Synonyms:	Proepiregulin Epiregulin (EPR) EREG
Observed MW:	19 kDa
Uniprot ID :	O14944

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab10012(EREG Antibody) at dilution of 1:100



A549 cells were subjected to SDS PAGE followed by western blot with FNab10012(EREG antibody) at dilution of 1:1000