

## PTGER4 antibody

### Product Information

Catalog No.:	FNab10018
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	PGE2 receptor EP4 subtype
Synonyms:	Prostaglandin E2 receptor EP4 subtype (PGE receptor EP4 subtype, PGE2 receptor EP4 subtype) Prostanoid EP4 receptor PTGER4 PTGER2
Observed MW:	50 kDa, 22 kDa
Uniprot ID :	P35408

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:2000
Image:	



Raji cells were subjected to SDS PAGE followed by western blot with FNab10018( PE2R4 Antibody) at dilution of 1:1000