

GK antibody

Product Information

Catalog No.:	FNab10019
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	Glycerol kinase
Synonyms:	Glycerol kinase (Glycerokinase) ATP:glycerol 3-phosphotransferase GK
Observed MW:	60 kDa
Uniprot ID :	P32189

Application

Reactivity:	Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:2000
Image:	



rat kidney were subjected to SDS PAGE followed by western blot with FNab10019(GLPK Antibody) at dilution of 1:1000