

APOBEC3A antibody

Product Information

Catalog No.:	FNab10020
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	APOBEC-3A
Synonyms:	DNA dC->dU-editing enzyme APOBEC-3A (A3A) Phorbolin-1 APOBEC3A
Observed MW:	26 kDa
Uniprot ID :	P31941

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab10020(ABC3A Antibody) at dilution of 1:200



rat testis were subjected to SDS PAGE followed by western blot with FNab10020(ABC3A Antibody) at dilution of 1:1000