

LMO1 antibody

Product Information

Catalog No.:	FNab10021
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	Rhombotin-1
Synonyms:	Rhombotin-1 Cysteine-rich protein TTG-1 LIM domain only protein 1 (LMO-1) T-cell translocation protein 1 LMO1 RBTN1 RHOM1 TTG1
Observed MW:	18 kDa
Uniprot ID :	P25800

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab10021(RBTN1 Antibody) at dilution of 1:200



rat brain were subjected to SDS PAGE followed by western blot with FNab10021(RBTN1 Antibody) at dilution of 1:1000