

LMO2 antibody

Product Information

Catalog No.:	FNab10022
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	Rhombotin-2
Synonyms:	Rhombotin-2 Cysteine-rich protein TTG-2 LIM domain only protein 2 (LMO-2) T-cell translocation protein 2 LMO2 RBTN2 RBTN1 RHOM2 TTG2
Observed MW:	17-25 kDa
Uniprot ID :	P25791

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded rat brain tissue slide using FNab10022(RBTN2 Antibody) at dilution of 1:200



k562 cells were subjected to SDS PAGE followed by western blot with FNab10022(RBTN2 Antibody) at dilution of 1:1000