

anti- AURKA antibody

Product Information

Catalog No.:	FNab10043
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	aurora kinase A
Synonyms:	AIK, ARK 1, ARK1, AURA, AURKA, Aurora A, aurora kinase A, Aurora related kinase 1, Aurora/IPL1 related kinase 1, AURORA2, Breast tumor amplified kinase, BTAK, hARK1, STK15, STK6, STK7
Observed MW:	46 kDa
Uniprot ID :	O14965

Application


Reactivity:	Human, Rat			
Tested Application:	ELISA, WB, IHC			
Recommended dilution: WB: 1:500 - 1:2000; IHC: 1:50 - 1:200				
Image:				

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com





Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab10043(AURKA Antibody) at dilution of 1:200

Hela cells were subjected to SDS PAGE followed by western blot with FNab10043(AURKA Antibody) at dilution of 1:1000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com