

TARDBP antibody

Product Information

Catalog No.:	FNab10050
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	TAR DNA binding protein
Synonyms:	TAR DNA-binding protein 43 (TDP-43) TARDBP TDP43
Observed MW:	43 kDa
Uniprot ID :	Q13148

Application

Reactivity:	Human,Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded rat kidney tissue slide using FNab10050(TADBP Antibody) at dilution of 1:200



rat liver were subjected to SDS PAGE followed by western blot with FNab10050(TADBP antibody) at dilution of 1:1000.