

PMP22 antibody

Product Information

Catalog No.:	FNab10061
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	peripheral myelin protein 22
Synonyms:	Peripheral myelin protein 22 (PMP-22) Growth arrest-specific protein 3 (GAS-3) PMP22 GAS3
Observed MW:	22 kDa
Uniprot ID :	Q01453

Application

Reactivity:	Human,Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab10061(PMP22 antibody) at dilution of 1:100



PC12 cells were subjected to SDS PAGE followed by western blot with FNab10061(PMP22 antibody) at dilution of 1:1000.