

## PPL antibody

### Product Information

|               |                                                                                                                 |
|---------------|-----------------------------------------------------------------------------------------------------------------|
| Catalog No.:  | FNab10072                                                                                                       |
| Size:         | 100µg                                                                                                           |
| Form:         | liquid                                                                                                          |
| Purification: | Immunogen affinity purified                                                                                     |
| Purity:       | ≥95% as determined by SDS-PAGE                                                                                  |
| Host:         | Rabbit                                                                                                          |
| Clonality:    | polyclonal                                                                                                      |
| Clone ID:     | None                                                                                                            |
| IsoType:      | IgG                                                                                                             |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Background


None

### Immunogen information

|              |                                                                                                               |
|--------------|---------------------------------------------------------------------------------------------------------------|
| Immunogen:   | periplakin                                                                                                    |
| Synonyms:    | Periplakin 190 kDa paraneoplastic pemphigus antigen 195 kDa cornified envelope precursor protein PPL KIAA0568 |
| Observed MW: | 204 kDa                                                                                                       |
| Uniprot ID : | O60437                                                                                                        |

### Application

|                       |                 |
|-----------------------|-----------------|
| Reactivity:           | Human           |
| Tested Application:   | ELISA, IHC      |
| Recommended dilution: | IHC: 1:20-1:200 |
| Image:                |                 |



Immunohistochemistry of paraffin-embedded human lung using FNab10072( Periplakin Antibody) at dilution of 1:100