

## PRKN antibody

### Product Information

Catalog No.:	FNab10088
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

None


### Immunogen information

Immunogen:	Parkinson disease(autosomal recessive, juvenile) 2, parkin
Synonyms:	E3 ubiquitin-protein ligase parkin (Parkin) Parkin RBR E3 ubiquitin-protein ligase Parkinson juvenile disease protein 2 (Parkinson disease protein 2) PRKN PARK2
Observed MW:	52 kDa
Uniprot ID :	O60260


### Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	





Immunohistochemistry of paraffin-embedded rat testis using FNab10088(Parkin antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab10088(Parkin antibody) at dilution of 1:1000