

ADCYAP1 antibody

Product Information

Catalog No.:	FNab10094
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	Pituitary adenylate cyclase-activating polypeptide
Synonyms:	Pituitary adenylate cyclase-activating polypeptide (PACAP) PACAP-related peptide Alternative names: PRP-48 Pituitary adenylate cyclase-activating polypeptide 27 (PACAP-27, PACAP27) Pituitary adenylate cyclase-activating polypeptide 38 (PACAP-38, PACAP38) ADCYAP1
Observed MW:	None
Uniprot ID :	P18509

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded rat brain tissue slide using FNab10094(PACA Antibody) at dilution of 1:200