

MED14 antibody

Product Information

Catalog No.:	FNab10101
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information


Immunogen:	Mediator of RNA polymerase II transcription subunit 14
Synonyms:	Mediator of RNA polymerase II transcription subunit 14 Activator-recruited cofactor 150 kDa component (ARC150) Cofactor required for Sp1 transcriptional activation subunit 2 (CRSP complex subunit 2) Mediator complex subunit 14 RGR1 homolog (hRGR1) Thyroid hormone receptor-associated protein complex 170 kDa component (Trap170) Transcriptional coactivator CRSP150 Vitamin D3 receptor-interacting protein complex 150 kDa component (DRIP150) MED14 ARC150 CRSP2 CXorf4 DRIP150 EXLM1 RGR1 TRAP170
Observed MW:	None
Uniprot ID :	O60244

Application

Reactivity:	Human
Tested Application:	ELISA, IHC

Recommended dilution: IHC: 1:50 - 1:500

Image:



Immunohistochemistry of paraffin-embedded human colon cancer using FNab10101(MED14 antibody) at dilution of 1:200