

## IL4R antibody

### Product Information

Catalog No.:	FNab10102
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	interleukin 4 receptor
Synonyms:	Interleukin-4 receptor subunit alpha (IL-4 receptor subunit alpha, IL-4R subunit alpha, IL-4R-alpha, IL-4RA) Soluble interleukin-4 receptor subunit alpha (Soluble IL-4 receptor subunit alpha, Soluble IL-4R-alpha, sIL4Ralpha/prot) Alternative names: IL-4-binding protein (IL4-BP) IL4R IL4RA
Observed MW:	35 kDa, 120 kDa
Uniprot ID :	P24394

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab10102(IL4RA antibody) at dilution of 1:200



K562 cells were subjected to SDS PAGE followed by western blot with FNab10102(IL4RA antibody) at dilution of 1:1000