

## MYF5 antibody

### Product Information

Catalog No.:	FNab10106
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	Myogenic factor 5
Synonyms:	Myogenic factor 5 (Myf-5) Class C basic helix-loop-helix protein 2 (bHLHc2) MYF5 BHLHC2
Observed MW:	None
Uniprot ID :	P13349

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50 - 1:500
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab10106(MYF5 antibody) at dilution of 1:200