

APOE antibody

Product Information

Catalog No.:	FNab10127
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	apolipoprotein E 33B
Synonyms:	Apolipoprotein E (Apo-E) APOE
Observed MW:	35 kDa
Uniprot ID :	P02649

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:50-1:200; IF: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human stomach tissue slide using FNab10127 (APOE Antibody) at dilution of 1:100.



SKOV3 cells were subjected to SDS PAGE followed by western blot with FNab10127 (APOE antibody) at dilution of 1:1000