

GRM2 antibody

Product Information

Catalog No.:	FNab10138
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	glutamate receptor, metabotropic 2
Synonyms:	Metabotropic glutamate receptor 2 (mGluR2) GRM2 GPRC1B MGLUR2
Observed MW:	96 kDa
Uniprot ID :	Q14416

Application

Reactivity:	Rat, Mouse
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNaB10138(GRM2 Antibody) at dilution of 1:200