

EWSR1 antibody

Product Information

Catalog No.:	FNab10142
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	Ewing sarcoma breakpoint region 1
Synonyms:	RNA-binding protein EWS EWS oncogene Ewing sarcoma breakpoint region 1 protein EWSR1 EWS
Observed MW:	37 kDa/68 kDa/100 kDa
Uniprot ID :	Q01844

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:50-1:200
Image:	



Immunofluorescent analysis of U2OS cells using FNab10142 (EWSR1 antibody) at dilution of 1:100



mouse spleen were subjected to SDS PAGE followed by western blot with FNab10142(EWSR1 antibody) at dilution of 1:1000