

GSDMB antibody

Product Information

Catalog No.:	FNab10147
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


The N-terminal moiety promotes pyroptosis. May be acting by homooligomerizing within the membrane and forming pores(PubMed:27281216). The physiological relevance of this observation is unknown(Probable).

Immunogen information

Immunogen:	gasdermin B
Synonyms:	Gasdermin-B Gasdermin-like protein Gasdermin-B, N-terminal (GSDMB-NT, p30) Gasdermin-B, C-terminal (GSDMB-CT, p16) GSDMB GSDML
Observed MW:	47 kDa
Uniprot ID :	Q8TAX9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution:	WB: 1:200-1:2000
Image:	



Hela cells were subjected to SDS PAGE followed by western blot with FNab10147(GSDMB Antibody) at dilution of 1:1000