

## TAB2 antibody

### Product Information

Catalog No.:	FNab10157
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Adapter linking MAP3K7/TAK1 and TRAF6. Promotes MAP3K7 activation in the IL1 signaling pathway. The binding of 'Lys-63'-linked polyubiquitin chains to TAB2 promotes autophosphorylation of MAP3K7 at 'Thr-187'. Involved in heart development.

### Immunogen information


Immunogen:	mitogen-activated protein kinase kinase kinase 7 interacting protein 2
Synonyms:	TGF-beta-activated kinase 1 and MAP3K7-binding protein 2 Mitogen-activated protein kinase kinase kinase 7-interacting protein 2 TAK1-binding protein 2 (TAB-2) TGF-beta-activated kinase 1-binding protein 2 TAB2 KIAA0733 MAP3K7IP2
Observed MW:	77 kDa
Uniprot ID :	Q9NYJ8

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded rat brain using FNab10157(TAB2 antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab10157(TAB2 Antibody) at dilution of 1:1000