

MAP2K7 antibody

Product Information

Catalog No.:	FNab10158
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	mitogen-activated protein kinase kinase 7
Synonyms:	Dual specificity mitogen-activated protein kinase kinase 7 (MAP kinase kinase 7, MAPKK 7) JNK-activating kinase 2 MAPK/ERK kinase 7 (MEK 7) Stress-activated protein kinase kinase 4 (SAPK kinase 4, SAPKK-4, SAPKK4) c-Jun N-terminal kinase kinase 2 (JNK kinase 2, JNKK 2) MAP2K7 JNKK2 MEK7 MKK7 PRKMK7 SKK4
Observed MW:	47 kDa, 51 kDa
Uniprot ID :	O14733

Application

Reactivity:	Human, rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded lung cancer using FNab10158(MKK7 antibody) at dilution of 1:100



Hela cells were subjected to SDS PAGE followed by western blot with FNab10158(MKK7 antibody) at dilution of 1:1000