

IRAK3 antibody

Product Information

Catalog No.:	FNab10159
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

This gene encodes a member of the interleukin-1 receptor-associated kinase protein family. Members of this family are essential components of the Toll/IL-R immune signal transduction pathways. This protein is primarily expressed in monocytes and macrophages and functions as a negative regulator of Toll-like receptor signaling. Mutations in this gene are associated with a susceptibility to asthma. Alternate splicing results in multiple transcript variants.


Immunogen information

Immunogen:	IRAK3
Synonyms:	Interleukin-1 receptor-associated kinase 3 (IRAK-3) IL-1 receptor-associated kinase M (IRAK-M) Inactive IL-1 receptor-associated kinase 3 IRAK3
Observed MW:	68 kDa
Uniprot ID :	Q9Y616

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000

Image:



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab10159(IRAK-3 antibody) at dilution of 1:1000