

C8A antibody

Product Information

Catalog No.:	FNab10168
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	Complement component C8 alpha chain
Synonyms:	Complement component C8 alpha chain Complement component 8 subunit alpha C8A
Observed MW:	65 kDa
Uniprot ID :	P07357

Application

Reactivity:	Human, rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded Rat liver tissue slide using FNab10168(CO8A antibody) at dilution of 1:200