

## TEAD1 antibody

### Product Information

Catalog No.:	FNab10169
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	Transcriptional enhancer factor TEF-1
Synonyms:	Transcriptional enhancer factor TEF-1 NTEF-1 Protein GT-IIC TEA domain family member 1 (TEAD-1) Transcription factor 13 (TCF-13) TEAD1 TCF13 TEF1
Observed MW:	50 kDa
Uniprot ID :	P28347

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:500
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab10169 (TEAD1 antibody) at dilution of 1:200



Various lysates were subjected to SDS PAGE followed by western blot with FNab10169 (TEAD1 antibody) at dilution of 1:1000