

CPE antibody

Product Information

Catalog No.:	FNab10172
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	Carboxypeptidase E
Synonyms:	Carboxypeptidase E (CPE) Carboxypeptidase H (CPH) Enkephalin convertase Prohormone-processing carboxypeptidase CPE CBPE
Observed MW:	53 kDa
Uniprot ID :	P16870

Application

Reactivity:	Human, rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded Rat brain tissue slide using FNab10172(CBPE antibody) at dilution of 1:200