

EZH2 antibody

Product Information

Catalog No.:	FNab10174
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	Histone-lysine N-methyltransferase EZH2
Synonyms:	Histone-lysine N-methyltransferase EZH2 ENX-1 Enhancer of zeste homolog 2 Lysine N-methyltransferase 6 EZH2 KMT6
Observed MW:	86 kDa
Uniprot ID :	Q15910

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab10174(EZH2 antibody) at dilution of 1:200