

## TSPAN12 antibody

### Product Information

Catalog No.:	FNab10176
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	Tetraspanin-12
Synonyms:	Tetraspanin-12 (Tspan-12) Tetraspan NET-2 Transmembrane 4 superfamily member 12 TSPAN12 NET2 TM4SF12
Observed MW:	22 kDa, 30-35 kDa
Uniprot ID :	O95859

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab10176(TSN12 antibody) at dilution of 1:200



Rat testis tissue were subjected to SDS PAGE followed by western blot with FNab10176(TSN12 antibody) at dilution of 1:1000