

TNC antibody

Product Information

Catalog No.:	FNab10178
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	Tenascin
Synonyms:	Tenascin (TN) Cytotactin GMEM GP 150-225 Glioma-associated-extracellular matrix antigen Hexabrachion JI Myotendinous antigen Neuronectin Tenascin-C (TN-C) TNC HXB
Observed MW:	241 kDa
Uniprot ID :	P24821

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab10178(TENA antibody) at dilution of 1:200