

SMARCB1 antibody

Product Information

Catalog No.:	FNab10186
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily B member 1
Synonyms:	SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily B member 1 BRG1-associated factor 47 (BAF47) Integrase interactor 1 protein SNF5 homolog (hSNF5) SMARCB1 BAF47 INI1 SNF5L1
Observed MW:	44 kDa
Uniprot ID :	Q12824

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab10186(SNF5 antibody) at dilution of 1:200