

## RPS6KA1 antibody

### Product Information

Catalog No.:	FNab10188
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	Ribosomal protein S6 kinase alpha-1
Synonyms:	Ribosomal protein S6 kinase alpha-1 (S6K-alpha-1) 90 kDa ribosomal protein S6 kinase 1 (p90-RSK 1, p90RSK1, p90S6K) MAP kinase-activated protein kinase 1a (MAPK-activated protein kinase 1a, MAPKAP kinase 1a, MAPKAPK-1a) Ribosomal S6 kinase 1 (RSK-1) RPS6KA1 MAPKAPK1A RSK1
Observed MW:	72-90 kDa
Uniprot ID :	Q15418

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab10188(KS6A1 antibody) at dilution of 1:200



HeLa cells were subjected to SDS PAGE followed by western blot with FNab10188(KS6A1 antibody) at dilution of 1:1000