

## EFS antibody

### Product Information

Catalog No.:	FNab10189
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	Embryonal Fyn-associated substrate
Synonyms:	Embryonal Fyn-associated substrate (hEFS) Cas scaffolding protein family member 3 EFS CASS3
Observed MW:	59 kDa
Uniprot ID :	O43281

### Application

Reactivity:	Mouse
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded mouse liver tissue slide using FNab10189(EFS antibody) at dilution of 1:200