

## DDR1 antibody

### Product Information

Catalog No.:	FNab10190
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	Epithelial discoidin domain-containing receptor 1
Synonyms:	Epithelial discoidin domain-containing receptor 1 (Epithelial discoidin domain receptor 1) CD167 antigen-like family member A Cell adhesion kinase Discoidin receptor tyrosine kinase HGK2 Mammary carcinoma kinase 10 (MCK-10) Protein-tyrosine kinase 3A Protein-tyrosine kinase RTK-6 TRK E Tyrosine kinase DDR Tyrosine-protein kinase CAK DDR1 CAK EDDR1 NEP NTRK4 PTK3A RTK6 TRKE
Observed MW:	97 kDa
Uniprot ID :	Q08345

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab10190(DDR1 antibody) at dilution of 1:200