

ANGPTL8 antibody

Product Information

Catalog No.:	FNab10195
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	Angiopoietin-like protein 8
Synonyms:	Angiopoietin-like protein 8 Betatrophin Lipasin Refeeding-induced fat and liver protein ANGPTL8 C19orf80 RIFL
Observed MW:	22 kDa
Uniprot ID :	Q6UXH0

Application

Reactivity:	Mouse
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded mouse liver tissue slide using FNab10195(ANGL8 antibody) at dilution of 1:200