

## FGA antibody

### Product Information

Catalog No.:	FNab10198
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	Fibrinogen
Synonyms:	Fibrinogen alpha chain Fibrinopeptide A Fibrinogen alpha chain FGA
Observed MW:	94 kDa
Uniprot ID :	P02671

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded mouse liver tissue slide using FNab10198(FIBA antibody) at dilution of 1:200