

## SHBG antibody

### Product Information

Catalog No.:	FNab10201
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	Sex hormone-binding globulin
Synonyms:	Sex hormone-binding globulin (SHBG) Sex steroid-binding protein (SBP) Testis-specific androgen-binding protein (ABP) Testosterone-estradiol-binding globulin (TeBG) Testosterone-estrogen-binding globulin SHBG
Observed MW:	32 kDa, 42 kDa, 44 kDa
Uniprot ID :	P04278

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney cancer tissue slide using FNab10201(SHBG antibody) at dilution of 1:200



Rat testis were subjected to SDS PAGE followed by western blot with FNab10201(SHBG antibody) at dilution of 1:1000