

BACE2 antibody

Product Information

Catalog No.:	FNab10202
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	Beta-secretase 2
Synonyms:	Beta-secretase 2 Aspartic-like protease 56 kDa Aspartyl protease 1 (ASP1, Asp 1) Beta-site amyloid precursor protein cleaving enzyme 2 (Beta-site APP cleaving enzyme 2) Down region aspartic protease (DRAP) Memapsin-1 Membrane-associated aspartic protease 1 Theta-secretase BACE2 AEPLC ALP56 ASP21
Observed MW:	56 kDa
Uniprot ID :	Q9Y5Z0

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab10202(BACE2 antibody) at dilution of 1:200