

## TOLLIP antibody

### Product Information

Catalog No.:	FNab10220
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	Toll-interacting protein
Synonyms:	Toll-interacting protein TOLLIP
Observed MW:	30 kDa
Uniprot ID :	Q9H0E2

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500; IF: 1:50-1:100
Image:	



Immunohistochemistry of paraffin-embedded mouse liver tissue slide using FNab10220(TOLLIP antibody) at dilution of 1:200



Immunofluorescent analysis of HepG2 cells using TOLLIP antibody(FNab10220) at 1:50 dilution.



Mcf-7 cells were subjected to SDS PAGE followed by western blot with FNab10220(TOLLIP antibody) at dilution of 1:1000