

AKT1 antibody

Product Information

Catalog No.:	FNab10227
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	RAC-alpha serine/threonine-protein kinase
Synonyms:	RAC-alpha serine/threonine-protein kinase Protein kinase B (PKB) Protein kinase B alpha (PKB alpha) Proto-oncogene c-Akt RAC-PK-alpha AKT1 PKB RAC
Observed MW:	48/55 kDa
Uniprot ID :	P31749

Application

Reactivity:	Human, rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded Rat kidney tissue slide using FNab10227(akt1 antibody) at dilution of 1:200



Hela cells were subjected to SDS PAGE followed by western blot with FNab10227(akt1 antibody) at dilution of 1:1000