

IKBKE antibody

Product Information

Catalog No.:	FNab10232
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	Inhibitor of nuclear factor kappa-B kinase subunit epsilon
Synonyms:	Inhibitor of nuclear factor kappa-B kinase subunit epsilon (I-kappa-B kinase epsilon, IKK-E, IKK-epsilon, IkbKE) Inducible I kappa-B kinase (IKK-i) IKBKE IKKE IKKI KIAA0151
Observed MW:	80 kDa
Uniprot ID :	Q14164

Application


Reactivity:	Human, rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500; IF: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded Rat kidney tissue slide using FNab10232(ikke antibody) at dilution of 1:200



Immunofluorescent analysis of U2OS cells using ikke antibody(FNab10232) at 1:50 dilution.



HL60 cells were subjected to SDS PAGE followed by western blot with FNab10232(ikke antibody) at dilution of 1:1000