
KRAS antibody

Product Information

Catalog No.:	FNab10239
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	GTPase KRas
Synonyms:	GTPase KRas K-Ras 2 Ki-Ras c-K-ras c-Ki-ras GTPase KRas, N-terminally processed KRAS KRAS2 RASK2
Observed MW:	21 kDa
Uniprot ID :	P01116

Application

Reactivity:	Human, rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200; IF: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded Rat kidney tissue slide using FNab10239(RASK antibody) at dilution of 1:100



A549 cells were subjected to SDS PAGE followed by western blot with FNab10239(rask antibody) at dilution of 1:1000