

RNF216 antibody

Product Information


Catalog No.:	FNab10259
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	E3 ubiquitin-protein ligase RNF216
Synonyms:	E3 ubiquitin-protein ligase RNF216 RING finger protein 216 RING-type E3 ubiquitin transferase RNF216 Triad domain-containing protein 3 Ubiquitin-conjugating enzyme 7-interacting protein 1 Zinc finger protein inhibiting NF-kappa-B RNF216 TRIAD3 UBCE7IP1 ZIN
Observed MW:	120-130 kDa
Uniprot ID :	Q9NWF9

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded rat kidney slide using FNab10259 (TRIAD3A Antibody) at dilution of 1:50



U87MG cells were subjected to SDS PAGE followed by western blot with FNab10259 (TRIAD3A Antibody) at dilution of 1:1000