

IRF7 antibody

Product Information


Catalog No.:	FNab10268
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	interferon regulatory factor 7
Synonyms:	Interferon regulatory factor 7 (IRF-7) IRF7
Observed MW:	56 kDa, 70 kDa
Uniprot ID :	Q92985

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab10268 (IRF7 Antibody) at dilution of 1:100



Mouse brain were subjected to SDS PAGE followed by western blot with FNab10268 (IRF7 Antibody) at dilution of 1:1000