

## ABCB11 antibody

### Product Information


Catalog No.:	FNab10278
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	ATP-binding cassette, sub-family B (MDR/TAP), member 11
Synonyms:	Bile salt export pump ATP-binding cassette sub-family B member 11 ABCB11 BSEP
Observed MW:	None
Uniprot ID :	O95342

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer using FNab10278(BSEP antibody) at dilution of 1:100